



Recombinant Methanocaldococcus jannaschii UPF0290 protein MJ1600 (MJ1600)

Product Code	CSB-CF701502MRU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q58995
Product Type	Transmembrane Protein
Immunogen Species	Methanocaldococcus jannaschii (strain ATCC 43067 / DSM 2661 / JAL-1 / JCM 10045 / NBRC 100440) (Methanococcus jannaschii)
Sequence	MFYRLLFASLWYILPAYVANASACIFGGGTPVDLGKNFIDGRRLLIGNGVTYRGC IFGILC GTLVGLIQGILVDFNIFNSLDFYGTVLDHVILAFFLSVGAIVGDAVGSGFIKRRRLNIE RGK PAPLLDQLDFVIGALAFGYIVAPIPYEMIIICLFTVFVHLLGNIIAYKLGKIDVWW
Source	in vitro E.coli expression system
Target Names	carS
Protein Names	Recommended name: UPF0290 protein MJ1600
Expression Region	1-177
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.