



Recombinant *Xenopus tropicalis* Transmembrane protein 196 (tmem196)

Product Code	CSB-CF023788XBF
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A4IIU3
Product Type	Transmembrane Protein
Immunogen Species	<i>Xenopus tropicalis</i> (Western clawed frog) (<i>Silurana tropicalis</i>)
Sequence	MCTSRQIIGSLLVLSVLEIGLGLSSVAVGAVTFFRVRMEQKPQLGDSSPVWSG ACFLICG VCGIFCAKKKSGLIMILFSACCICGLIGGILNIQFLRALNKRSSSLYSLYLASMSLA CIG ISGCTISTWLTCRLASYEQRRMFLEREHSLHHSHEMTEKDTENITNGGGPLAL NGRV
Source	in vitro E.coli expression system
Target Names	tmem196
Protein Names	Recommended name: Transmembrane protein 196
Expression Region	1-177
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.