



Recombinant Sorghum bicolor CASP-like protein Sb05g025790 (Sb05g025790)

Product Code	CSB-CF513653FJS
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	C5Y7C6
Product Type	Transmembrane Protein
Immunogen Species	Sorghum bicolor (Sorghum) (Sorghum vulgare)
Sequence	MAEEVWKALSLLFRIAALGLSLAAAIVMATASQLVIGGGGGHESSESSYSVSFGQ YNALRYF VAAGAISAVCSAAALYLFAVRADFTVVVVSLPLVPVLDAAAQGLFSAAGAAFA TRDVVG GGTSAGRGSSVCDAAGAFCGRVTVAAAVCAFAAVSVATAALASRDAGGGSS EGRRFEW
Source	in vitro E.coli expression system
Target Names	Sb05g025790
Protein Names	Recommended name: CASP-like protein Sb05g025790
Expression Region	1-178
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.