



Recombinant Arabidopsis thaliana CASP-like protein At2g39530 (At2g39530)

Product Code	CSB-CF811476DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8GWD5
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MAPPPPAPPSVTLRTVLLLLRVLTA AFLITVVLISTNTVTLEISSTS IKLPFNDVY AYR YMLSAAVIGLVYAVVQLFLTISQFATGKTHPLTYQDFYGDKVISYLLATGSAAG FGVSK DLKDTYIALIEFDSTDPVDKFFSKGYASASLLLFAFVSLAVLSVFSSLALS KR PVP VS
Source	in vitro E.coli expression system
Target Names	At2g39530
Protein Names	Recommended name: CASP-like protein At2g39530
Expression Region	1-178
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.