



Recombinant Escherichia fergusonii Probable intracellular septation protein A (yciB)

Product Code	CSB-CF484197EOR
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7LS24
Product Type	Transmembrane Protein
Immunogen Species	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Sequence	MKQFLDFLPLVVFFAFYKIYDIYAATAALIVATAIVLIYSWVRFRKVEKMALITFVL VVV FGGLTLFFHNDEFIKWKVTVIYALFAGALLVSQWVMKKPLIQRMLGKELSLPQP VWSKLN LAWAVFFILCGLANIYIAFWLPQNIWVNFKVFGTLTALTLIFTLLSGVYIYRHLPQD DKS
Source	in vitro E.coli expression system
Target Names	yciB
Protein Names	Recommended name: Probable intracellular septation protein A
Expression Region	1-179
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.