



Recombinant *Gossypium hirsutum* CASP-like protein XL3 (XL3)

Product Code	CSB-CF692296GHB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5K4H9
Product Type	Transmembrane Protein
Immunogen Species	<i>Gossypium hirsutum</i> (Upland cotton) (<i>Gossypium mexicanum</i>)
Sequence	MELSIQKIEALIRLSTIVMLVLTACLIGLDSQTKVIFYVQKKASFKDLRALVGLLYIT SL AAAYNLLQLCCSSFSASYKGTSLQSYAYLAWLRYILDQAVVYAVFAGNLALE HSFLVLT GEENFQWLKWCNKYTRFCTQIGGSLLCGFVASLLMFSIASISAFNLFRQYSPT KFMHLKL
Source	in vitro E.coli expression system
Target Names	XL3
Protein Names	Recommended name: CASP-like protein XL3
Expression Region	1-180
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.