



Recombinant Mpv17-like protein (T18D3.9)

Product Code	CSB-CF800087CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q7YWV6
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MVILFIRRLATNPLSTQMCIAGTISGSGDCLAQYLSHNQEWDRWRTARFSFL SSCFMA PSLFIWFRLLLEKVKGNKSLLLVKKLCIDQLCFSPCFNAAILFNLRLQLHQSAEK SWDLL KEDWFNIYATSLKVWPFVQVVNLCFVPLNYRVILNQVVAFFWNCYLSYITQKPI DHIEQF Y
Source	in vitro E.coli expression system
Target Names	T18D3.9
Protein Names	Recommended name: Mpv17-like protein
Expression Region	1-181
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.