



Recombinant Human Plasmolipin (PLLIP)

Product Code	CSB-CF896501HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9Y342
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAEFPSKVSTRTSSPAQGAEASVSALRPDLGFVRSRLGALMLLQLVLGLLVWA LIADTPY HLYPAYGWVMFVAVFLWLVTIVLFLNLYLFQLHMKLYMVPWPLVLMIFNISATVL YITAFI ACSAAVDLTSLRGTRPYNQRAAASFFACLVMIAYGVSAFFSYQAWRGVGSNA ATSQMAGG YA
Source	in vitro E.coli expression system
Target Names	PLLIP
Protein Names	Recommended name: Plasmolipin Alternative name(s): Plasma membrane proteolipid
Expression Region	1-182
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.