



Recombinant *Oryza sativa* subsp. japonica CASP-like protein Os11g0649700 (Os11g0649700, LOC_Os11g42970)

Product Code	CSB-CF650865OFG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q2R0D1
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Sequence	MSYGCQVSDDEPNGSKAVSLLLRLSTLALALTSAVVMATASECTVVQLNGVVA TITYKDF PPFVYLVGFNIAAAMLEAAAIYLRSTGGGDDDDDEGFKGKLPGILLVVIDVAVQA LVYTA TGGAFAAVSAYGPQINACGAGAGRFCGQVHQSKLLSFAGSAAVGLAVVFRDV SLPFSLWP TSSD
Source	in vitro E.coli expression system
Target Names	Os11g0649700
Protein Names	Recommended name: CASP-like protein Os11g0649700
Expression Region	1-184
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.