



Recombinant *Oryza sativa* subsp. japonica Casparian strip membrane protein Os06g0231050 (Os06g0231050, LOC_Os06g12500)

Product Code	CSB-CF727840OFG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q67X40
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Sequence	MEGSEEHGETSKAPLSRGVSKGVSILDVILRFVAIIGTLASAIAMGTTNQTLPFF TQFIR FKAQYSDLPTLTFFVANSIVSAYLILSLPLSIVHVIRSRKYSRLILIFFDAAMLAL VT AGASAAAIVYLAHKGNARANWLAICQQFDSFCERISGSLIGSFAAMVVLVLLIF LSAIA LARR
Source	in vitro E.coli expression system
Target Names	Os06g0231050
Protein Names	Recommended name: Casparian strip membrane protein Os06g0231050
Expression Region	1-184
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.