



Recombinant Zea mays CASP-like protein 10

Product Code	CSB-CF469621ZAX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B6SZA7
Product Type	Transmembrane Protein
Immunogen Species	Zea mays (Maize)
Sequence	MAAAARVSEVKAEGLLRGACAALAAAAALLVGLSTQTETVLLVRKKATVKDVQ ALWVLAM AAAAAAGYHLLQLLKCLYLGRVGGARPCRRSSRALAWTCLLLDKACAYTTFAT TVAAAQA CVVALDGAHAVQWTKLCNIYTRFCEQVAGSLVLGMLAAVGTAVLSAASARNV FRHYSSLE TYAAH
Source	in vitro E.coli expression system
Protein Names	Recommended name: CASP-like protein 10
Expression Region	1-185
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.