



Recombinant Pongo abelii Translocon-associated protein subunit gamma (SSR3)

Product Code	CSB-CF719255PYX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5RCD7
Product Type	Transmembrane Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Sequence	MAPKGSCKQQSEEDLLLQDFSRNLSAKSSALFFGNFIVSAIPIWLYWRIWHM DLIQSAV LYSVMTLVSTYLVAFAFKNVKFLKHKVAQKREDAVSKEVTRKLSEADNRKMS RKEKDER ILWKKNEVADYEATTFSIFYNNTLFLVVIVASFFILKNFNPTVNYILSISASSGLIA LL STGSK
Source	in vitro E.coli expression system
Target Names	SSR3
Protein Names	Recommended name: Translocon-associated protein subunit gamma Short name= TRAP-gamma Alternative name(s): Signal sequence receptor subunit gamma Short name= SSR-gamma
Expression Region	1-185
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.