



Recombinant Human Prenylated Rab acceptor protein 1 (RABAC1)

Product Code	CSB-CF887054HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9UI14
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAAQKDQQKDAEAEGLSGTLLPKLIPSGAGREWLERRRATIRPWSTFVDQQ RFSRPRNL GELCQRLVRNVEYYQSNEYVFLGLILYCVVTSPMLLVALAVFFGACYILYLRTL ESKLV LFGREVSPAHQYALAGGISFPFFWLAGAGSAVFWVLGATLVVIGSHAAFHQIE AVDGEEL QMEPV
Source	in vitro E.coli expression system
Target Names	RABAC1
Protein Names	Recommended name: Prenylated Rab acceptor protein 1 Alternative name(s): PRA1 family protein 1
Expression Region	1-185
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.