



Recombinant Zea mays CASP-like protein 12

Product Code	CSB-CF478390ZAX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B6TUW9
Product Type	Transmembrane Protein
Immunogen Species	Zea mays (Maize)
Sequence	MGAAGLKVPEMALRLCVVPLSLASLWEMASNAQADDTYGEVKFSDLSGFSYL VGVNAVTA AYAVASVLASSFKRPLAARYDWVVLVMDQASAYLLVTSASAAAELLQLARHGD RGVSWGE ACSYFGRFCGKATVSLALHAAALACFAALSLVSAFRVFSSRCHPPPDADGQPP KHARDEE QRVYHY
Source	in vitro E.coli expression system
Protein Names	Recommended name: CASP-like protein 12 Alternative name(s): Protein salicylic acid-induced 1
Expression Region	1-186
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.