



Recombinant Escherichia fergusonii UPF0059 membrane protein yebN (yebN)

Product Code	CSB-CF483291EOR
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7LPM5
Product Type	Transmembrane Protein
Immunogen Species	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Sequence	MNITATVLLAFGMSMDAFAASIGKGATLHKPKFSEALRTGLIFGAVETLTPLIGW GMGML ASRFVLEWNHWIAFVLLIFLGGRMIIEGIRGGDEDEEPRRRHGFWLLVTTAIA TSLDAM AVGVGLAFLQVNIATALAIGCATLIMSTLGMMVGRFIGPILGKKAIEILGGLVLIGI GVQ ILWTHFHG
Source	in vitro E.coli expression system
Target Names	mntP
Protein Names	Recommended name: UPF0059 membrane protein yebN
Expression Region	1-188
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.