



Recombinant Escherichia coli UPF0059 membrane protein yebN (yebN)

Product Code	CSB-CF631136EGW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q1RAW9
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain UTI89 / UPEC)
Sequence	MNITATVLLAFGMSMDAFAASIGKGATLHKPKFSEALRTGLIFGAVETLTPLIGW GMGML ASRFVLEWNHWIAFVLLIFLGGRMIIEGFRGADDEDEEPRRRHGFWLLVTTAIA TSLDAM AVGVGLAFLQVNIATALAIGCATLIMSTLGMMVGRFIGSIIGKKAIEILGGLVLIGI GVQ ILWTHFHG
Source	in vitro E.coli expression system
Target Names	mntP
Protein Names	Recommended name: UPF0059 membrane protein yebN
Expression Region	1-188
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.