



Recombinant Danio rerio Transmembrane protein 205 (tmem205)

Product Code	CSB-CF023799DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A1L2F6
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MATEGDPTDFVKVLHLLVISFTWGMQVWVSFIAGFVLISQVSMHTFGLVQSKL FPVYFYC LLGNAVSLAVYAVYHPRELLDWHEGIQMSLFFVAVIMAGLNAQWFGPSATE NMLVMQEI EKEHGLGNQVGMSSNREGYTKLREQDPKYKEHRSTFYRYHGLSNLCNLIGFF CITVNLIIY LALNLGTI
Source	in vitro E.coli expression system
Target Names	tmem205
Protein Names	Recommended name: Transmembrane protein 205
Expression Region	1-188
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.