



Recombinant Danio rerio Transmembrane protein 229b (tmem229b)

Product Code	CSB-CF600568DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q08CG9
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MATTVTPEPLTALSRWYLYAIHGYFCEVMFTAAWEFVVNCNWKFPVTSVWA LFIYGTCI LIVERMYLCLKDRCNVLLRCIIYTLWTFWFEFGTGFLLRQFNACPWDYSEFKYN FMGLIT AEYAVPWFCASFIVERLVIRNTLRLRFDEVAESGQAEERLDRGGGGRRGGRRG RGARAGAT SANGYVKVD
Source	in vitro E.coli expression system
Target Names	tmem229b
Protein Names	Recommended name: Transmembrane protein 229b
Expression Region	1-189
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.