



Recombinant *Saccharomyces cerevisiae* Uncharacterized membrane protein YGR016W (YGR016W)

Product Code	CSB-CF344863SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P53209
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MSRLRRFNRKILSLSSDYTHDGESDQEDVSILPLDTEEQEELIQKFETNAHITNK LYINL LSILYLLYGGLLMILVRKSRGYIKLALLAGANSLICSCITLRYDIVNDYLLFKKFKL RVS NFSINIINIILLVLMAWISFNHVVEDKKTVLCLQVPMFLFWVAVLVKRWARNIEDE IADL RCLKYKYKNA
Source	in vitro E.coli expression system
Target Names	YGR016W
Protein Names	Recommended name: Uncharacterized membrane protein YGR016W
Expression Region	1-190
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.