



Recombinant Arabidopsis thaliana CASP-like protein At4g15630 (At4g15630)

Product Code	CSB-CF809255DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8L8Z1
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MEHESKNKVDGMEMEKGGKESGSRKGLELTMRVLALVLTMVAATVLGVAKQT KVVPIKLI PTLPPLNVSTTAKASYLSAFVYNISANAIACGYTAISIVIVMISKGKRKSKLLMAVL IGD LMMVALLFSSTGAAGAIGLMGRHGNKHMWKKVCGVFGKFCNQAAVSVAILI ASVVFML LVVLDALKLP
Source	in vitro E.coli expression system
Target Names	At4g15630
Protein Names	Recommended name: CASP-like protein At4g15630
Expression Region	1-190
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.