



Recombinant *Selaginella moellendorffii* CASP-like protein SELMODRAFT_75865 (SELMODRAFT_75865)

Product Code	CSB-CF520435SYF
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	D8QNV6
Product Type	Transmembrane Protein
Immunogen Species	<i>Selaginella moellendorffii</i> (Spikemoss)
Sequence	MGAYDGAEAPRAAPASTAANSRPSRLLLLHSLLLRLVAVVVSILVIAVMVHAKQ RVMIFK AEWDNSKAFVALVAISAICLGYSFLQFILSAFHLCSSKSWKSPTKCWAWMNFIA DQILTYA MLGAAAAAAELAYIAKNGSSRAQWQPICSTFNTFCTRAGASIILSFI AVLALANS SAISA YHLFRRPSSSV
Source	in vitro E.coli expression system
Target Names	SELMODRAFT_75865
Protein Names	Recommended name: CASP-like protein SELMODRAFT_75865
Expression Region	1-191
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.