



Recombinant *Oryza sativa* subsp. *japonica*

Probable protein-S-isoprenylcysteine O-methyltransferase (ICMT)

Product Code	CSB-CF773725OFG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q7XSR9
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Sequence	MAARAQAWLFAAALVIFHGSEYVLAAAFHGRRNVTATSLISKQYVLAMSFAM LEHLTEA LLFPELKEYWVFSYVGLVMVIIGEVIKRLAVVTAGRSFTHVIRIHVEDQHKLITHG VYRL MRHPGYSGFLIWA VGTQVMLCNPLSTVAFTLVLRFFSKRIPYEEFFLRQFFG REYEEYA QKVHSGLPFIE
Source	in vitro E.coli expression system
Target Names	ICMT
Protein Names	Recommended name: Probable protein-S-isoprenylcysteine O-methyltransferase EC= 2.1.1.100 Alternative name(s): Isoprenylcysteine carboxymethyltransferase Prenylated protein carboxyl methyltransferase Prenylcysteine carboxyl methyltr
Expression Region	1-191
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.