



Recombinant Glycine max CASP-like protein 7

Product Code	CSB-CF513097GGV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	C6SVQ5
Product Type	Transmembrane Protein
Immunogen Species	Glycine max (Soybean) (Glycine hispida)
Sequence	MASTDKPGGDPEYRTSSTPAPAGVDYFKFDVILRFLFFAASLVAVVVIVTANQT EVIRVP QPVPWPAKFRYSPAFVYFVAALSVTGLYSIITTLASLLASNKPALKTKLLLYFILW DALI LGIIASATGTAGGVAYLGLKGNRHVVGWNKICHVYDKFCRHVGASIAVALFGS VVTVLLI WLSAYSIHSRVPK
Source	in vitro E.coli expression system
Protein Names	Recommended name: CASP-like protein 7
Expression Region	1-193
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.