



Recombinant Arabidopsis thaliana CASP-like protein At5g54980 (At5g54980)

Product Code	CSB-CF872212DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9FFT2
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MRDNNNNNTREEERSSSSKQQQPQAPMSLKIIDSCLRLSVVPLSVATIWLTVT NHESNPD YGNLEYNSIMGLKYMVGVS AISAIYALLSTVSSWVTCLVSKAWLFFIPDQVLAY VMTTSV AGATEIVYLLNKGDKIVTWSEMCSYPHYCSKLTIALGLHVFVLFVFLFLSVISAY RAFS PFDPPCDSQTNNDA
Source	in vitro E.coli expression system
Target Names	At5g54980
Protein Names	Recommended name: CASP-like protein At5g54980
Expression Region	1-194
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.