



Recombinant Arabidopsis thaliana Protein RER1B (RER1B)

Product Code	CSB-CF527392DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O48671
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MEGSGGDSGSMATPVQKKVHEAWRVYQYYLDKTTPHSTNRWIGTLVFFLIYC LRVYSIHG FYIISYGLGIYLLNLLIGFLSPLVDPELEVSDGATLPTRGSDEFKPFIRRLPEFKF WYSM TKAFCIAFLMTFFSVFDVPVFWPILLCYWVVLVLTMRRIAAMIKHKYIPFSIGK QKYS GRKSSANSGGGSRAD
Source	in vitro E.coli expression system
Target Names	RER1B
Protein Names	Recommended name: Protein RER1B Short name= AtRER1B
Expression Region	1-195
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.