



# Recombinant Human Mpv17-like protein (MPV17L)

|                          |                                                                                                                                                                                                                                                                                    |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-CF650065HU                                                                                                                                                                                                                                                                     |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.                                                                                                                                                                                                                  |
| <b>Uniprot No.</b>       | Q2QL34                                                                                                                                                                                                                                                                             |
| <b>Product Type</b>      | Transmembrane Protein                                                                                                                                                                                                                                                              |
| <b>Immunogen Species</b> | Homo sapiens (Human)                                                                                                                                                                                                                                                               |
| <b>Sequence</b>          | MAGWWPALSRAARRHPWPTNVLLYGSLVSAGDALQQRLLQGREANWRQTRR<br>VATLVVTFHA<br>NFNYYVWLRLLEALPGRAPHALLAKLLCDQVVGAPIAVSAFYVGMSSILQGKDDI<br>FLDLKQ<br>KFWNTYLSGLMYWPFVQLTNFSLVPVQWRTAYAGVCGFLWATFICFSQQSG<br>DGTFKSAFT ILYTKGTSATEGYPKK                                                        |
| <b>Source</b>            | in vitro E.coli expression system                                                                                                                                                                                                                                                  |
| <b>Target Names</b>      | MPV17L                                                                                                                                                                                                                                                                             |
| <b>Protein Names</b>     | Recommended name: Mpv17-like protein Alternative name(s): M-LP homolog<br>Short name= M-LPH                                                                                                                                                                                        |
| <b>Expression Region</b> | 1-196                                                                                                                                                                                                                                                                              |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                      |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |