



Recombinant Arabidopsis thaliana Vacuolar iron transporter homolog 2 (At1g76800)

Product Code	CSB-CF890410DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9SRD3
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MDQSGSNTNMDIEKESTTFDYSKRSQWLRAAVLGANDGLVSTASLMMGVGA VKHDVKAMI LSGFAGMVAGACSM AIGEFVSVYSQYDIEVAQMERDSVEIEKEKLPSPMQAAA ASALAFS AGAIVPLLAAAFVKEYKMRIISVVVAVTV ALMVFGWLGAALGKAPAVRSSARVL FGWLA MAVTFGLTKLIGLYGL
Source	in vitro E.coli expression system
Target Names	At1g76800
Protein Names	Recommended name: Vacuolar iron transporter homolog 2 Alternative name(s): Protein NODULIN-LIKE 2
Expression Region	1-196
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.