



Recombinant *Saccharomyces cerevisiae* Autophagy-related protein 33 (ATG33)

Product Code	CSB-CF453637SVP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B3RHM5
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain RM11-1a) (Baker's yeast)
Sequence	MSVCLAITKGIAVSSIGLYSGLLASASLITSTTPLEVLGTGSLTPTLTTLKNAATALG AFA STFFCVSFFGAPPSLRHPYLLYGMLAAPLSSFVLGCASNYQSRKYSKVSKESS LPEDSK PAASELSDSIIDLGEDNHASENTPRDGKPAATTVSKPAEALHTGPPIHTKNLIAA TAIAI VGFVQAVIGVYGEGQFI
Source	in vitro E.coli expression system
Target Names	ATG33
Protein Names	Recommended name: Autophagy-related protein 33
Expression Region	1-197
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.