



Recombinant Arabidopsis thaliana Protein-S-isoprenylcysteine O-methyltransferase A (ICMTA)

Product Code	CSB-CF866990DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9FMW9
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MTEIFS DTSIRQLSQMLLSLIFFHISEYILAITIHGASNVTLSSLLITKHYALAMLLS LL EYLTEIILFPGLKQHWWVSNFGLIMIIVGEIIRKAAITAGRSFTHLIKINYEEHHGL VT HGVYRLMRHPSYCGFLIWSVGTQVMLCNPVSAVAFVWVRFFAQRIPIYEEY FLNQFFGV QYLEYAESVASGVPFVN
Source	in vitro E.coli expression system
Target Names	ICMTA
Protein Names	Recommended name: Protein-S-isoprenylcysteine O-methyltransferase A Short name= AtICMTA EC= 2.1.1.100 Alternative name(s): Isoprenylcysteine carboxylmethyltransferase A Prenylated protein carboxyl methyltransferase A Prenylcy
Expression Region	1-197
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.