



Recombinant *Oryza sativa* subsp. japonica CASP-like protein Os03g0817100 (Os03g0817100, LOC_Os03g60250)

Product Code	CSB-CF497861OFG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B9F6Z0
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Sequence	MSSSGPPAGDGRDDASGPGPAGAAAAADGSVPVSRIVERWKMEPAAARAR LLRAVAWL FLLALVVMASNKHGHGGAQDFDNYPEYTYCLGISIIAVLYTTAQVTRDVHRLS WGRDVI AGRKAAAVVDFAGDQVVAYLLMSALSAAAPVTDYMRQAADNLFTDSAAAAS MAFLAFLA AGLSALVSGYNLAMEVLV
Source	in vitro E.coli expression system
Target Names	Os03g0817100
Protein Names	Recommended name: CASP-like protein Os03g0817100
Expression Region	1-198
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.