



Recombinant *Oryza sativa* subsp. *japonica* CASP-like protein Os12g0514300 (Os12g0514300, LOC_Os12g32970)

Product Code	CSB-CF610259OFG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q0IN16
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Sequence	MAAAMGLERKAKVAEVALRCAVCAALAAALVGTGSQTRTFFSLEKKARFTD MKALVLL VAAHGAAVYSLLQLARCAAAA AWKGG SNGGA AVVAWSVFSCDQAVAYALM AATAAALQS SVVGKRGQPELQWMPVCGLYGAFRRVGEGLAAAVAAGLAAVLLAAVSAFNL FRLYGGGG GGRKSSAGAVSGNGANTW
Source	in vitro <i>E.coli</i> expression system
Target Names	Os12g0514300
Protein Names	Recommended name: CASP-like protein Os12g0514300
Expression Region	1-198
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.