



Recombinant Arabidopsis thaliana Vacuolar iron transporter homolog 4 (At3g43660)

Product Code	CSB-CF888900DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9M2C0
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MESNNNNLNLDMEKDQETTFDYSKRAQWLRAAVLGANDGLVSTASLMMGIG AVKQDVRIM LLTGFAGLVAGACSM AIGEFISVYSQYDIEVAQMKRESGGETKKEKLPSTQA AIASALA FTLGAIVPLLA AAFVKEYKVRIGVIVA AAVTLALVMFGWLGAVLGKAPVVKSLVRV LIGGW LAMAITFGFTKLVGSHGL
Source	in vitro E.coli expression system
Target Names	At3g43660
Protein Names	Recommended name: Vacuolar iron transporter homolog 4 Alternative name(s): Protein NODULIN-LIKE 4
Expression Region	1-198
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.