



Recombinant Escherichia coli Putative undecaprenyl-diphosphatase ybjG (ybjG)

Product Code	CSB-CF301157ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P75806
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	MLENLNLSLFSLINATPDSAPWMISLAIFIADLITVPLLAVVLWLWGLTAQRQL VIKI AIALAVSLFVSWTMGHLFPHDRPFVENIGYNFLHHAADDSFPSDHGTVIFTFAL AFLCWH RLWSGSLLMVLAVVIAWSRVYLGVHWPLDMLGGLLAGMIGCLSAQIIWQAMG HKLYQRLQ SWYRVCFALPIRKGWVRD
Source	in vitro E.coli expression system
Target Names	ybjG
Protein Names	Recommended name: Putative undecaprenyl-diphosphatase ybjG EC= 3.6.1.27 Alternative name(s): Undecaprenyl pyrophosphate phosphatase
Expression Region	1-198
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.