



Recombinant Escherichia fergusonii Sulfoxide reductase heme-binding subunit YedZ (yedZ)

Product Code	CSB-CF486747EOR
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7LRM9
Product Type	Transmembrane Protein
Immunogen Species	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Sequence	MRLTAKQVIWLKLVCLHLAGFLPFVWLAWAINHGGLSADPVKDIQHFTGRTALK FLLATLL ITPLARYAKQPLLIRTRRLLGLWCFAWATLHLTSYALLELGVNNLALFGQELITR PYLTL GIVSWLILLALAFTSTQAMQRNLGKHWQQLHNFVYLVAILAPIHYLWSVKIISPQ PILYA GLALLLLALRYKKFRTWLK
Source	in vitro E.coli expression system
Target Names	msrQ
Protein Names	Recommended name: Sulfoxide reductase heme-binding subunit YedZ Alternative name(s): Flavocytochrome YedZ
Expression Region	1-199
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.