



Recombinant *Saccharomyces cerevisiae* Putative uncharacterized protein YIR017W-A (YIR017W-A)

Product Code	CSB-CF312744SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q03885
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MFRLMSFNFKDLISLMSLSFAFQNSFSLFNCSIRSRSLICVFKFCSLFMFSK FFCFL RMRNRCDASVFFLLRSCLILSSSLEVLTGASSSFTTTAVVAATFPFLGSRSLN CETSWP LFIMLMSSSALPLLTSSNDSKREPSSSESLDIRFCRCSESLFTFIPFLAMSMFVD FFSHP CFALILTNTICLRNYHWL
Source	in vitro E.coli expression system
Target Names	YIR017W-A
Protein Names	Recommended name: Putative uncharacterized protein YIR017W-A
Expression Region	1-199
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.