



Recombinant *Oryza sativa* subsp. *indica* CASP-like protein OsI_06397 (OsI_06397)

Product Code	CSB-CF386065OFF
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A2X2I0
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Sequence	MRSGEGSTAAAAAAEEEEKVKVAAPFRLAELGLRVCAVPLAVASVWEMATNKQ VDETYGEV RFSDLSGFRYLVWINAITAAYSVASILLSSCRFITRFDWLIFILDQASAYLLLTSAS AAA EVVYLAREGDREVSWSWGEVCSYFGRFCGAATVSVALNAAALLCFMALSLSIFR VFTKFNP PSQSNSKQQLSQEQGKPVVSG
Source	in vitro E.coli expression system
Target Names	OsI_06397
Protein Names	Recommended name: CASP-like protein OsI_06397
Expression Region	1-201
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.