



Recombinant *Oryza sativa* subsp. *indica* CASP-like protein Osl_07795 (Osl_07795)

Product Code	CSB-CF489613OFF
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B8AEC9
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Sequence	MEAGEEIEDGEPSTPTYKAHHPPPHLPPPMRSSGVSLVLSVADLVLRVFAIGG TAGSAIA MATTSETLPFAAPFVRFRAEYSDLPTLMFFVVAASSVVCAYLVLSLPASVVHVVR PGARSS RAILAFLDTVMLALLTASASAAAAIVYLAHRGSARANWLGICQQFTSFCQRITAS LVGSF AAVVLVALVFLSALSARRA
Source	in vitro E.coli expression system
Target Names	Osl_07795
Protein Names	Recommended name: CASP-like protein Osl_07795
Expression Region	1-201
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.