



Recombinant *Oryza sativa* subsp. *japonica* CASP-like protein Os02g0219900 (Os02g0219900, LOC_Os02g12760)

Product Code	CSB-CF744363OFG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6YW53
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Sequence	MRSGEGSTAAAAAAEEEEKVKVAAPFRLAELGLRVCAVPLAVASVWEMATNKQ VDETYGEV RFSDLSGFRYLVWINAITAAYSVASILLSSCRFITRFDWLIFLLDQASAYLLLTSA SAAA EVVYLAREGDREVSWGEVCSYFGRFCGAATVSVALNAAALLCFMALSLISAFR VFTKFNP PSQSNSKQQLSQEQGKPVVSG
Source	in vitro <i>E.coli</i> expression system
Target Names	Os02g0219900
Protein Names	Recommended name: CASP-like protein Os02g0219900
Expression Region	1-201
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.