



Recombinant Human Transmembrane 4 L6 family member 1 (TM4SF1)

Product Code	CSB-CF023615HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P30408
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MCYGKCARCIGHSLVGLALLCIAANILLYFPNGETKYASENHLSRFVWFFSGIV GGLLM LLPAFVFIGLEQDDCCGCCGHENCGKRCAMLSSVLAALIGIAGSGYCVIVAALG LAEGPL CLDSLQWNYTFASTEGQYLLDTSTWSECTEPKHIVEWNVSLFSILLALGGIEF ILCLIQ VINGVLGGICGFCCSHQQQYDC
Source	in vitro E.coli expression system
Target Names	TM4SF1
Protein Names	Recommended name: Transmembrane 4 L6 family member 1 Alternative name(s): Membrane component chromosome 3 surface marker 1 Tumor-associated antigen L6
Expression Region	1-202
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.