



Recombinant Selaginella moellendorffii CASP-like protein SELMODRAFT_272229 (SELMODRAFT_272229)

Product Code	CSB-CF518700SYF
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	D8T829
Product Type	Transmembrane Protein
Immunogen Species	Selaginella moellendorffii (Spikemoss)
Sequence	MSEHRIPVAADKQISPPISAGEQKGCKGLKRTDMLRFAAFVCCAVTMVVLITD KQTSAI QVPGFNLTITKTVSFDLAKAFVYLVSAAGIGAGYTLLVLVLSIISAERSKAIAWFI FVF DQLITYVLLAAAAASTEVAVMGAHAPPEASWLKVC SLFGRFCHQLGASLVTSFI STVLFA FSAAISAYYLF SNTNVRPAYS KG
Source	in vitro E.coli expression system
Target Names	SELMODRAFT_272229
Protein Names	Recommended name: CASP-like protein SELMODRAFT_272229
Expression Region	1-203
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.