



Recombinant Escherichia coli Cytochrome o ubiquinol oxidase subunit 3 (cyoC)

Product Code	CSB-CF364658ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0ABJ3
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	MATDTLTHATAHAHEHGHHDAGGTKIFGFWIYLMSDCILFSILFATYAVLVNGT AGGPTG KDIFELPFVLVETFLLLFSSITYGMAAIAMYKNNKSQVISWLALTWLFGAGFIGM EIYEF HHLIVNGMGPDRSGFLSAFFALVGTHGLHVTSGLIWMAVLMVQIARRGLTSTN RTRIMCL SLFWHFLDVVWICVFTVVYLMGAM
Source	in vitro E.coli expression system
Target Names	cyoC
Protein Names	Recommended name: Cytochrome o ubiquinol oxidase subunit 3 EC= 1.10.3.- Alternative name(s): Cytochrome o ubiquinol oxidase subunit III Ubiquinol oxidase chain C
Expression Region	1-204
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.