



Recombinant *Oryza sativa* subsp. *japonica* CASP-like protein Os05g0344400 (Os05g0344400, LOC_Os05g27790)

Product Code	CSB-CF722802OFG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5W6M3
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Sequence	MSAAVAASSGAPAADVEKGAAAADANVDGGGAPAAAAASGEGVVS AVVRRW RRQDLLEKS GSALRVAAWAFSLLAFVVMGANDHGDWRQFEHYEEYRYVVAIGVLAFIYTTLQ LVRHGVR LTGGQDLQGKVAVLVDFAGDQVTAYLLMSAVSAAIPITNRMREGADNVFTDSS AASISMA FFAFLCLALSALVSGFKLAKQTYI
Source	in vitro <i>E.coli</i> expression system
Target Names	Os05g0344400
Protein Names	Recommended name: CASP-like protein Os05g0344400
Expression Region	1-204
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.