



Recombinant Human Putative claudin-24 (CLDN24)

Product Code	CSB-CF413094HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A6NM45
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MALIFRTAMQSVGLLLSLLGWILSIITTYLPHWKNLNLDLNEMENWTMGLWQTC VIQEEV GMQCKDFDSFLALPAELRVSRILMFLSNGLGFLGLLVSGFGLDCLRIGESQRDL KRRLLI LGGILSWASGITALVPVSWVAHKTQVEFDENVPDFVPRWEFGEALFLGWFA GLSLLLGG CLLNCAACSSHAPLALGHYAVAQMQ
Source	in vitro E.coli expression system
Target Names	CLDN24
Protein Names	Recommended name: Putative claudin-24
Expression Region	1-205
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.