



# Recombinant Human Transmembrane protein 150B (TMEM150B)

|                          |                                                                                                                                                                                                                                                                                    |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-CF405636HU                                                                                                                                                                                                                                                                     |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.                                                                                                                                                                                                                  |
| <b>Uniprot No.</b>       | A6NC51                                                                                                                                                                                                                                                                             |
| <b>Product Type</b>      | Transmembrane Protein                                                                                                                                                                                                                                                              |
| <b>Immunogen Species</b> | Homo sapiens (Human)                                                                                                                                                                                                                                                               |
| <b>Sequence</b>          | VTNRTVDLSKGFYISICGSFPPQSCIFSQVLNMGAALAAWICIVRYHQLRDWG<br>VRRWPN<br>QLILWTGLLCALGTSVVGNFQEKNQRPTHLAGAFLAFILGNVYFWLQLLLWRL<br>KRLPQPG<br>AAWIGPLRLGLCSVCTILIVAMIVLHACSLRSVSAACEWVVAMLLFALFGLLAVD<br>FSALE SCTL CVQPWPSLSPPPASPISLPVQL                                                |
| <b>Source</b>            | in vitro E.coli expression system                                                                                                                                                                                                                                                  |
| <b>Target Names</b>      | TMEM150B                                                                                                                                                                                                                                                                           |
| <b>Protein Names</b>     | Recommended name: Transmembrane protein 150B Alternative name(s):<br>Transmembrane protein 224                                                                                                                                                                                     |
| <b>Expression Region</b> | 28-233                                                                                                                                                                                                                                                                             |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                      |
| <b>Protein Length</b>    | Full Length of Mature Protein                                                                                                                                                                                                                                                      |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |