



Recombinant Danio rerio Transmembrane protein 150B (tmem150b)

Product Code	CSB-CF656297DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q32PK2
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	IAVSNNSVNITIEFPYISTCGAYTPQSCLFAQICNICCVLALWIVVIRFQQIRDLGR SSH LNTAGLVLGFISSIGISILGNFQQTIIQEVHLLGALMAFFLGLAYFWIQAFITYFSP PSR DNKWLVPVRFVLC SQCTCMVICMFVLHSTGFRSAAAICEWILVMCFFALFGVF AAEFRHI DFHKLTVQKEGLKVANNDNVVWTVQDVQ
Source	in vitro E.coli expression system
Target Names	tmem150b
Protein Names	Recommended name: Transmembrane protein 150B Alternative name(s): Transmembrane protein 224
Expression Region	25-232
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.