



Recombinant Macaca mulatta (Rhesus macaque) Platelet-activating factor receptor (PTAFR)

Product Code	CSB-CF018946MOW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P35366
Product Type	Transmembrane Protein
Immunogen Species	Macaca mulatta (Rhesus macaque)
Sequence	MEPHDSSHVDSEFRYTLFPPIVYSIIFVLGVIANGYVLWVWFARLYPSKKFNEIKIFM VNLT MADMLFLITLPLWIVYYQNGGNWIFPKFLCNLAGCLFFINTYCSVAFLGVITYNR FQAVT RPIKTAQANTRKRGISLSLVIWVAIVGAASYFFILDSTNTVPNSAGSGNITRCFE HYEKG SVPVLIHIFIVFSFFLVFLIILFCNLV
Source	in vitro E.coli expression system
Target Names	PTAFR
Protein Names	Recommended name: Platelet-activating factor receptor Short name= PAF-R Short name= PAFr
Expression Region	1-208
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.