



Recombinant *Oryza sativa* subsp. *indica* CASP-like protein Osl_15195 (Osl_15195, OSIGBa0131F24.1)

Product Code	CSB-CF488328OFF
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B8ARW3
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Sequence	MSKMAEEKVLAAPATVDGGMQSSGDLQASSAAAARVRPVETLLRAAPLGLCV AAMAIMLR NSVTNEYGTVSYSDLGGFKYLVYANGLCAAYSLASAFYIAVPRPATLSRSWV FLLDQVF TYLILAAGAASAELLYLAYNGDKEVTWSEACGVFGGFCRQARTSVAITFASVA CYILLSL ISSYRLFSAYDPPQPQLGNKGVEIAAFPR
Source	in vitro <i>E.coli</i> expression system
Target Names	Osl_15195
Protein Names	Recommended name: CASP-like protein Osl_15195
Expression Region	1-209
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.