



Recombinant Vitis vinifera CASP-like protein GSVIVT00034332001 (VIT_09s0002g03780)

Product Code	CSB-CF416556VFQ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A7P756
Product Type	Transmembrane Protein
Immunogen Species	Vitis vinifera (Grape)
Sequence	MNGLKTPPEIGIQLPEAKVAAETGTMSGPLVPPRSRSDRSVRRGTDVAHVVLRFV CLLTSVI ALSLMATAKEAASISIYGFLLPVSSKWSFSDSFEYLVGVSAAVAHAHALLQLIISVS RLLR KSPVIPSRNHAWLIFAGDQAFAYAMLSAGSAASGVTNLNRTGIRHSPLPNFCK PLRSFCD HVAASIAFTFFSCFLLATSAILDVIWLSKY
Source	in vitro E.coli expression system
Target Names	VIT_09s0002g03780
Protein Names	Recommended name: CASP-like protein GSVIVT00034332001
Expression Region	1-210
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.