



Recombinant Dictyostelium discoideum Protein YIPF5 homolog (yipf5)

Product Code	CSB-CF709803DKK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q54QY3
Product Type	Transmembrane Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Sequence	MNNNSFNFIDSQYSTPQGAYYDNTGRMGGGGGMGGPTDSFDNELPLLEEL GINFDHIRS KTLNVLNPLKKIDSHIMDDTDLGGPILFGLLLGFSLLMSGKIQFGYIYGLGLIGCV SMYI VLNLMSEKGDIDYRVISVLGYCLLPMIFLSFTSLIININGMVGYILIGFAIWWSTYSA SK MFVKVLSMIDQRILVAYPVGLLYTGFAITAF
Source	in vitro E.coli expression system
Target Names	yipf5
Protein Names	Recommended name: Protein YIPF5 homolog
Expression Region	1-212
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.