



Recombinant Danio rerio Tetraspanin-31 (tspan31)

Product Code	CSB-CF767229DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q7ZUB3
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MVCGGFTCSKNALCSLNVVYMLVGLLLIVVAAWGKGFIVSSIHIGGVIAVGFF LQLIA IVGLIGAVHHHQVMLFFYMVILFVVFLLFQFGVSCSCLAMNQGGQEKLESSWKI MSNDTR ISLEKKLDCCGLFNSTNLQADIMSDLHLCTSPCTQKKECVTCGLKMLQHSSEAL KILGGV GLFFSFTEILGVWLAMRYRNQKDPRANPSAFL
Source	in vitro E.coli expression system
Target Names	tspan31
Protein Names	Recommended name: Tetraspanin-31 Short name= Tspan-31 Alternative name(s): Sarcoma-amplified sequence homolog
Expression Region	1-212
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.